



SygmaLift

### High performance skin lifting in combination therapy



ISO  
13485

Certification

CE 0029

CE Certification

- Pulsed H.F.U. Focal Ultrasounds, allowing treatment from first layer to deeper tissue.
- Intense beam of cold laser 635nm for skin tightening.
- Treatment performed in about 40min.
- Touch screen with smart technology philosophy.
- Affordable, no hidden extra costs.
- Ergonomic transducers.
- Data base & Special Leasing option.
- Compact and mobile unit
- Weight 38 kilo
- 3 handpieces
- 200/240 Volts - 50/60 Hz.



MedixSysteme

info@g3lasers.com

G3 Lasers  
empowering practices™

+91 98408 92502



SygmaLift

### Exclusive HIFU Lines Probe

Providing faster treatment, greater precision and a unique solution for sub mental fat, facial lifting and rejuvenation.

SygmaLift®

- Focal lines of US for modeling
- Micro focal ultrasounds for lifting
- Cold laser for toning and glow

Photos copyright : MedixSysteme - Y. Nakamura - M. Jung



EXCLUSIVE  
HIFU LINES



### High Intensity Focal Ultrasounds

High Fractional Focused ultrasound, to concentrate up to 1000 intersecting beams of ultrasound energy with extreme precision on a target depth in the body. Much like a magnifying glass can focus multiple beams of light on a single point. Where each individual beam of focused ultrasound passes through the tissue, it has a progressive thermal effect down to the adipose tissue to stimulate fat lyses.

### Circular Micro Focal Ultrasounds

Three Micro Focal US ultrasound, to stimulate in triple focal pulsed mode at various deepness the skin for regeneration of tissue. Focal Ultrasounds have been shown clinically to be an effective noninvasive method for the stimulation of bioactivity. Secondly H.F.U. rejuvenation leads to denaturation of the collagen by mechanically heating the dermis and stimulation of collagen synthesis. Results are progressive in time.

### Intensive Cold Laser

Scanning Diode Laser 635nm, Lasers are based upon a principle of physics, wavelengths. To be considered "therapeutic" they must have a wavelength between 630 and 905 nanometer (nm). The strength of the wavelength determines how deep into the tissue the laser can reach. Lower level lasers are used for superficial skin treatment, while higher wavelengths stimulate deeper cells and tissue.



### SyigmaLift Unique high tech facelift solutions Combination treatment that works

#### Performance on :

- Chin / Jaw lines reshaping
- Eye contour / Eye Brow lifting
- Skin Glow, Tightening and Rejuvenation
- Fine lines



#### Long lasting results

- Synthesis of fibroblasts
- Increased oxygenation of tissue
- Stimulation of collagen growth

#### Fast results with no downtime

- Specific area procedures
- Maintenance session after 6 months recommended
- Walk in walk out procedure lasting 40 minutes

#### Cost effective Solution

- No consumables, except some specific gels
- Gold ceramics for longevity hand pieces
- Internationally patented, Made in France

#### Ergonomic

- User friendly interface, guiding physicians for protocol
- Three handpieces, each conforming to specific application
- Hair stimulation protocol

